

## Heptanoic acid methyl ester

## Chemical Properties

CAS No. : 106-73-0


Formula: C<sub>8</sub>H<sub>16</sub>O<sub>2</sub>

Molecular Weight: 144.21

Pure form: -20°C for 3 years | In solvent: -80°C for 1

Storage: year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Heptanoic acid methyl ester is a fatty acid methyl ester widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.9343 mL	34.6717 mL	69.3433 mL
5 mM	1.3869 mL	6.9343 mL	13.8687 mL
10 mM	0.6934 mL	3.4672 mL	6.9343 mL
50 mM	0.1387 mL	0.6934 mL	1.3869 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481