


4-Ethynylpyridine

Chemical Properties

CAS No. :	2510-22-7
Formula:	C ₇ H ₅ N
Molecular Weight:	103.124
Storage:	Store at low temperature Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	4-Ethynylpyridine, with CAS No. 2510-22-7, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 4-Ethynylpyridine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Solubility Information

Solubility	DMSO: 20.62 mg/mL (199.95 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	9.6974 mL	48.4872 mL	96.9744 mL
5 mM	1.9395 mL	9.6974 mL	19.3949 mL
10 mM	0.9697 mL	4.8487 mL	9.6974 mL
50 mM	0.1939 mL	0.9697 mL	1.9395 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481