


## [3-[(4-Methylpiperidino)methyl]phenyl]methanamine

## Chemical Properties

CAS No. :	690632-06-5
Formula:	C14H22N2
Molecular Weight:	218.3379
Storage:	Pure form: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	[3-[(4-Methylpiperidino)methyl]phenyl]methanamine, with CAS No. 690632-06-5, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. [3-[(4-Methylpiperidino)methyl]phenyl]methanamine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.580 mL	22.9001 mL	45.8001 mL
5 mM	0.916 mL	4.580 mL	9.160 mL
10 mM	0.458 mL	2.290 mL	4.580 mL
50 mM	0.0916 mL	0.458 mL	0.916 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481