



## Solubility Information

Solubility	DMSO: 10 mg/mL (23.5 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 1 mg/mL (2.35 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3503 mL	11.7514 mL	23.5029 mL
5 mM	0.4701 mL	2.3503 mL	4.7006 mL
10 mM	0.235 mL	1.1751 mL	2.3503 mL
50 mM	0.047 mL	0.235 mL	0.4701 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

- Gellibert F, et al. J Med Chem. 2006, 49(7), 2210-2221.  
 Petersen M, et al. Kidney Int, 2008, 73(6), 705-715.  
 Tan SM, et al. Am J Physiol Heart Circ Physiol, 2010, 298(5), H1415-1425.  
 Lagares D, et al. Arthritis Rheum, 2010, 62(3), 878-889.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481