

Compound N025-0034

Chemical Properties


CAS No. :

Formula: C48H78O18

Molecular Weight: 943.134

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N025-0034 is a useful organic compound for research related to life sciences and the catalog number is T131312.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0603 mL	5.3015 mL	10.603 mL
5 mM	0.2121 mL	1.0603 mL	2.1206 mL
10 mM	0.106 mL	0.5301 mL	1.0603 mL
50 mM	0.0212 mL	0.106 mL	0.2121 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481