


Biotin-PEG6-amine

Chemical Properties

CAS No. :

Formula: C₂₄H₄₆N₄O₈S

Molecular Weight: 550.71



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Biotin-PEG6-amine is a PEG-based linker for PROTACs, joining two essential ligands crucial for forming PROTAC molecules and enabling selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands connected by a linker: one targets an E3 ubiquitin ligase, and the other binds the target protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8158 mL	9.0792 mL	18.1584 mL
5 mM	0.3632 mL	1.8158 mL	3.6317 mL
10 mM	0.1816 mL	0.9079 mL	1.8158 mL
50 mM	0.0363 mL	0.1816 mL	0.3632 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481