


Polygalin C

Chemical Properties

CAS No. :	934768-05-5
Formula:	C ₂₈ H ₃₂ O ₁₆
Molecular Weight:	624.54
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Polygalin C is a natural product for research related to life sciences. The catalog number is TN6529 and the CAS number is 934768-05-5.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6012 mL	8.0059 mL	16.0118 mL
5 mM	0.3202 mL	1.6012 mL	3.2024 mL
10 mM	0.1601 mL	0.8006 mL	1.6012 mL
50 mM	0.032 mL	0.1601 mL	0.3202 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481