

## 3-oxoandrost-4-en-17-yl hexanoate

## Chemical Properties


CAS No. :

Formula: C<sub>25</sub>H<sub>38</sub>O<sub>3</sub>

Molecular Weight: 386.576

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-oxoandrost-4-en-17-yl hexanoate is a useful organic compound for research related to life sciences and the catalog number is T131615.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5868 mL	12.9339 mL	25.8679 mL
5 mM	0.5174 mL	2.5868 mL	5.1736 mL
10 mM	0.2587 mL	1.2934 mL	2.5868 mL
50 mM	0.0517 mL	0.2587 mL	0.5174 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481