


BLXA-4

Chemical Properties

CAS No. : 362516-29-8
 Formula: C₂₃H₃₄O₅
 Molecular Weight: 390.51
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	BLXA-4 is used potentially to treat gingival diseases.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5608 mL	12.8038 mL	25.6075 mL
5 mM	0.5122 mL	2.5608 mL	5.1215 mL
10 mM	0.2561 mL	1.2804 mL	2.5608 mL
50 mM	0.0512 mL	0.2561 mL	0.5122 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Van Dyke TE, Hasturk H, Kantarci A, Freire MO, Nguyen D, Dalli J, Serhan CN. Proresolving nanomedicines activate bone regeneration in periodontitis. J Dent Res. 2015 Jan;94(1):148-56. doi: 10.1177/0022034514557331. PubMed PMID: 25389003; PubMed Central PMCID: PMC4270812.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481