


2',3',5'-Tri-O-benzoyl-3-deazauridine

Chemical Properties

CAS No. : 23220-74-8
 Formula: C₃₁H₂₅N₃O₉
 Molecular Weight: 555.53
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Nucleosides and Reagents - 3-Dezauridine; Pyridine nucleoside
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8001 mL	9.0004 mL	18.0008 mL
5 mM	0.360 mL	1.8001 mL	3.6002 mL
10 mM	0.180 mL	0.900 mL	1.8001 mL
50 mM	0.036 mL	0.180 mL	0.360 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481