


## Daturaturin A aglycone

## Chemical Properties

CAS No. :	904665-71-0
Formula:	C <sub>28</sub> H <sub>38</sub> O <sub>5</sub>
Molecular Weight:	454.607
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Daturaturin A aglycone is a natural product for research related to life sciences. The catalog number is TN3775 and the CAS number is 904665-71-0.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1997 mL	10.9984 mL	21.9969 mL
5 mM	0.4399 mL	2.1997 mL	4.3994 mL
10 mM	0.220 mL	1.0998 mL	2.1997 mL
50 mM	0.044 mL	0.220 mL	0.4399 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481