

2-{5-[2-({2-[5-(2-hydroxybutyl)oxolan-2-yl]propanoyl}oxy)butyl]oxolan-2-yl}propanoic acid

Chemical Properties

CAS No. :

Formula: C22H38O7

Molecular Weight: 414.54

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	2-{5-[2-({2-[5-(2-hydroxybutyl)oxolan-2-yl]propanoyl}oxy)butyl]oxolan-2-yl}propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4123 mL	12.0616 mL	24.1231 mL
5 mM	0.4825 mL	2.4123 mL	4.8246 mL
10 mM	0.2412 mL	1.2062 mL	2.4123 mL
50 mM	0.0482 mL	0.2412 mL	0.4825 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481