

Compound T130688

Chemical Properties

CAS No. :

Formula: C₂₆H₃₆O₆

Molecular Weight: 444.57

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	methyl (3aR,3bR,9aR,11aS)-2-hydroxy-2,3b,6,6,9a,11a-hexamethyl-11-methylidene-1,3,7-trioxo-hexadecahydro-1H-cyclopenta[a]phenanthrene-3a-carboxylate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2494 mL	11.2468 mL	22.4936 mL
5 mM	0.4499 mL	2.2494 mL	4.4987 mL
10 mM	0.2249 mL	1.1247 mL	2.2494 mL
50 mM	0.045 mL	0.2249 mL	0.4499 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481