

2-[2-(2,3-dihydroxy-3-methylbutyl)-4-hydroxyphenoxy]-6-(hydroxymethyl)oxane-3,4,5-triol

Chemical Properties


CAS No. :

Formula: C17H26O9

Molecular Weight: 374.39

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-[2-(2,3-dihydroxy-3-methylbutyl)-4-hydroxyphenoxy]-6-(hydroxymethyl)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.671 mL	13.3551 mL	26.7101 mL
5 mM	0.5342 mL	2.671 mL	5.342 mL
10 mM	0.2671 mL	1.3355 mL	2.671 mL
50 mM	0.0534 mL	0.2671 mL	0.5342 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481