

## (10Z,13Z)-15,16-dihydroxyoctadeca-10,13-dienoic acid

## Chemical Properties

CAS No. :

Formula: C<sub>18</sub>H<sub>32</sub>O<sub>4</sub>

Molecular Weight: 312.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(10Z,13Z)-15,16-dihydroxyoctadeca-10,13-dienoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2005 mL	16.0026 mL	32.0051 mL
5 mM	0.6401 mL	3.2005 mL	6.401 mL
10 mM	0.3201 mL	1.6003 mL	3.2005 mL
50 mM	0.064 mL	0.3201 mL	0.6401 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481