


2'-C-Methyl -2',3',5'-tri-O-benzoyluridine

Chemical Properties

CAS No. :	23643-36-9
Formula:	C ₃₁ H ₂₆ N ₂ O ₉
Molecular Weight:	570.55
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2'-C-Methyl -2',3',5'-tri-O-benzoyluridine is a 2'-C-Methyl nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7527 mL	8.7635 mL	17.5269 mL
5 mM	0.3505 mL	1.7527 mL	3.5054 mL
10 mM	0.1753 mL	0.8763 mL	1.7527 mL
50 mM	0.0351 mL	0.1753 mL	0.3505 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481