


## 1H-1,2,4-Triazole-3-carboxylic acid

## Chemical Properties

CAS No. :	4928-87-4
Formula:	C3H3N3O2
Molecular Weight:	113.07
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	It has been granted to have property to inhibit protein isoprenyl transferases.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 30 mg/mL (265.32 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.8441 mL	44.2204 mL	88.4408 mL
5 mM	1.7688 mL	8.8441 mL	17.6882 mL
10 mM	0.8844 mL	4.422 mL	8.8441 mL
50 mM	0.1769 mL	0.8844 mL	1.7688 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Lindner M, et al. Sci Pharm. 2010; 78(2):195-214.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481