

Farfugin A

Chemical Properties

CAS No. :	36061-18-4
Formula:	C ₁₅ H ₁₈ O
Molecular Weight:	214.30
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Farfugin A is a natural product that can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6664 mL	23.3318 mL	46.6636 mL
5 mM	0.9333 mL	4.6664 mL	9.3327 mL
10 mM	0.4666 mL	2.3332 mL	4.6664 mL
50 mM	0.0933 mL	0.4666 mL	0.9333 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481