


16-Oxolyclanitin-29-yl p-coumarate

Chemical Properties

CAS No. : 140701-70-8
 Formula: C₃₉H₅₄O₈
 Molecular Weight: 650.84
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description 16-Oxolyclanitin-29-yl p-coumarate is a natural product for research related to life sciences. The catalog number is TN5971 and the CAS number is 140701-70-8.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5365 mL	7.6824 mL	15.3648 mL
5 mM	0.3073 mL	1.5365 mL	3.073 mL
10 mM	0.1536 mL	0.7682 mL	1.5365 mL
50 mM	0.0307 mL	0.1536 mL	0.3073 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481