


D-Pyrohomoglutamic Acid

Chemical Properties

CAS No. :	72002-30-3
Formula:	C ₆ H ₉ NO ₃
Molecular Weight:	143.14
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	D-Pyrohomoglutamic acid is a building block. It has been used in the synthesis of immunoproteasome low molecular mass polypeptide 2 (LMP2) subunit inhibitors and is a component of substance P peptide analogs.
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.9862 mL	34.9308 mL	69.8617 mL
5 mM	1.3972 mL	6.9862 mL	13.9723 mL
10 mM	0.6986 mL	3.4931 mL	6.9862 mL
50 mM	0.1397 mL	0.6986 mL	1.3972 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481