

(2S)-1-O-p-coumaroyl-2-O-acetyl-3-O-β-D-glucopyranosylglycerol

Chemical Properties

CAS No. :

Formula: C₂₀H₂₆O₁₁

Molecular Weight: 442.42

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(2S)-1-O-p-coumaroyl-2-O-acetyl-3-O-β-D-glucopyranosylglycerol is a useful research compound for many research applications, including the fields of chemistry and life sciences, among others.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 50 mg/mL (113.01 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2603 mL	11.3015 mL	22.603 mL
5 mM	0.4521 mL	2.2603 mL	4.5206 mL
10 mM	0.226 mL	1.1301 mL	2.2603 mL
50 mM	0.0452 mL	0.226 mL	0.4521 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481