


## Dihydrosenkyunolide C

## Chemical Properties

CAS No. :	195142-72-4
Formula:	C <sub>12</sub> H <sub>14</sub> O <sub>3</sub>
Molecular Weight:	206.24
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Dihydrosenkyunolide Cis a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 100 mg/mL (484.87 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8487 mL	24.2436 mL	48.4872 mL
5 mM	0.9697 mL	4.8487 mL	9.6974 mL
10 mM	0.4849 mL	2.4244 mL	4.8487 mL
50 mM	0.097 mL	0.4849 mL	0.9697 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Qiaohua Qin, et al. Design, synthesis and biological evaluation of novel indolin-2-one derivatives as potent second-generation TRKs inhibitors. Eur J Med Chem. 2023 May 5;253:115291.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481