


ARD-353

## Chemical Properties

CAS No. : 561068-32-4  
 Formula: C32H39N3O4  
 Molecular Weight: 529.67  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	ARD-353 is a nonpeptide delta receptor agonist with cardioprotective effects.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.888 mL	9.4398 mL	18.8797 mL
5 mM	0.3776 mL	1.888 mL	3.7759 mL
10 mM	0.1888 mL	0.944 mL	1.888 mL
50 mM	0.0378 mL	0.1888 mL	0.3776 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481