

Compound T130603

Chemical Properties


CAS No. :

Formula: C44H68O19

Molecular Weight: 901.01

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

3a-hydroxy-7-[(4-methoxy-6-methyl-5-[[3,4,5-trihydroxy-6-((hydroxymethyl)oxan-2-yl]oxy)methyl)oxan-2-yl]oxy]oxan-2-yl]oxy]-9a,11a-dimethyl-1-(5-oxo-2,5-dihydrofuran-3-yl)-hexadecahydro-1H-cyclopenta[a]phenanthren-2-yl acetate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1099 mL	5.5493 mL	11.0987 mL
5 mM	0.222 mL	1.1099 mL	2.2197 mL
10 mM	0.111 mL	0.5549 mL	1.1099 mL
50 mM	0.0222 mL	0.111 mL	0.222 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481