


## Sugeroside

## Chemical Properties

CAS No. :	41743-57-1
Formula:	C <sub>26</sub> H <sub>42</sub> O <sub>8</sub>
Molecular Weight:	482.61
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sugeroside is a bitter-tasting diterpenoid glycoside, primarily isolated from the leaves of <i>Ilex sugerokii</i> var., which exhibits inhibitory activity against A549 cells (IC <sub>50</sub> = 2.39 μM).
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0721 mL	10.3603 mL	20.7207 mL
5 mM	0.4144 mL	2.0721 mL	4.1441 mL
10 mM	0.2072 mL	1.036 mL	2.0721 mL
50 mM	0.0414 mL	0.2072 mL	0.4144 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Kawakami S, et al. Eight ent-Kaurane Diterpenoid Glycosides Named Diosmariosides A-H from the Leaves of *Diospyros maritima* and Their Cytotoxic Activity. *Chem Pharm Bull (Tokyo)*. 2018;66(11):1057-1064.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481