

Anti-MERS-3A1 mAb

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Anti-MERS-3A1 mAb (MERS-3A1) is a high-affinity human IgG1 monoclonal antibody derived from CHO cells that effectively inhibits the interaction between the MERS-CoV spike protein and the DPP4 receptor[1].
Targets(IC50)	SARS-CoV

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481