


## 4',4''-Dihydroxylanigorootin

## Chemical Properties

CAS No. :	448949-11-9
Formula:	C <sub>38</sub> H <sub>22</sub> O <sub>8</sub>
Molecular Weight:	606.586
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4',4''-Dihydroxylanigorootin is a natural product from Musa cv."Thepanom" (BBB).
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6486 mL	8.2428 mL	16.4856 mL
5 mM	0.3297 mL	1.6486 mL	3.2971 mL
10 mM	0.1649 mL	0.8243 mL	1.6486 mL
50 mM	0.033 mL	0.1649 mL	0.3297 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481