


Wilforlide A acetate

Chemical Properties

CAS No. :	84104-80-3
Formula:	C ₃₂ H ₄₈ O ₄
Molecular Weight:	496.72
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Wilforlide A (T1) serves as a quality control standard of Tripterygium Glycosides that be listed in Drug Standard of Ministry of Public Health of the People's Republic of China. T1 has anti-inflammatory effects.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0132 mL	10.066 mL	20.1321 mL
5 mM	0.4026 mL	2.0132 mL	4.0264 mL
10 mM	0.2013 mL	1.0066 mL	2.0132 mL
50 mM	0.0403 mL	0.2013 mL	0.4026 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481