

Hydroxy-PEG1-C2-methyl ester

Chemical Properties

CAS No. : 93673-82-6


Formula: C₆H₁₂O₄

Molecular Weight: 148.16

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Hydroxy-PEG1-C2-methyl ester, a polyethylene glycol (PEG) derived linker, is widely used in the production of PROTACs (proteolysis targeting chimeras). [1]
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together: one binds to an E3 ubiquitin ligase, and the other targets a specific protein. By leveraging the intracellular ubiquitin-proteasome system, PROTACs facilitate the selective degradation of target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7495 mL	33.7473 mL	67.4946 mL
5 mM	1.3499 mL	6.7495 mL	13.4989 mL
10 mM	0.6749 mL	3.3747 mL	6.7495 mL
50 mM	0.135 mL	0.6749 mL	1.3499 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481