


FK-143

Chemical Properties

CAS No. : 163136-03-6
 Formula: C₄₀H₄₄N₂O₃
 Molecular Weight: 600.79
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	FK-143 is a 5-alpha reductase inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6645 mL	8.3224 mL	16.6448 mL
5 mM	0.3329 mL	1.6645 mL	3.329 mL
10 mM	0.1664 mL	0.8322 mL	1.6645 mL
50 mM	0.0333 mL	0.1664 mL	0.3329 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481