


3,5-Dicaffeoyl-epi-quinic acid

Chemical Properties

CAS No. :	879305-14-3
Formula:	C ₂₅ H ₂₄ O ₁₂
Molecular Weight:	516.45
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3,5-Dicaffeoyl-epi-quinic acid exhibits anti-inflammatory activity and may improve mast cell-mediated inflammatory diseases.
Targets(IC50)	Caspase,Histamine Receptor,IL Receptor,JNK,p38 MAPK,TNF

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9363 mL	9.6815 mL	19.363 mL
5 mM	0.3873 mL	1.9363 mL	3.8726 mL
10 mM	0.1936 mL	0.9681 mL	1.9363 mL
50 mM	0.0387 mL	0.1936 mL	0.3873 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481