


MDL-72392

Chemical Properties

CAS No. : 85278-68-8
 Formula: C₉H₁₀FNO
 Molecular Weight: 167.18
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	MDL-72392 is a monoamine oxidase inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9816 mL	29.9079 mL	59.8158 mL
5 mM	1.1963 mL	5.9816 mL	11.9632 mL
10 mM	0.5982 mL	2.9908 mL	5.9816 mL
50 mM	0.1196 mL	0.5982 mL	1.1963 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481