


## Songoroside A

## Chemical Properties

|                   |   |
|-------------------|---|
| CAS No. :         | 61617-29-6  |
| Formula:          | C <sub>35</sub> H <sub>56</sub> O <sub>7</sub>  |
| Molecular Weight: | 588.81  |
| Storage:          | Powder: -20°C for 3 years   In solvent: -80°C for 1 year<br>Actual storage temperature shall be subject to the COA. |



## Biological Description

|               |   |
|---------------|---|
| Description   | Songoroside A is a natural product of Sanguisorba, Rosaceae. The catalog number is T7801 and the CAS number is 61617-29-6. Songoroside A can be used as a reference standard. |
| Targets(IC50) | Others  |

## Preparing Stock Solutions

|       | 1mg       | 5mg       | 10mg       |
|-------|-----------|-----------|------------|
| 1 mM  | 1.6983 mL | 8.4917 mL | 16.9834 mL |
| 5 mM  | 0.3397 mL | 1.6983 mL | 3.3967 mL  |
| 10 mM | 0.1698 mL | 0.8492 mL | 1.6983 mL  |
| 50 mM | 0.034 mL  | 0.1698 mL | 0.3397 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481