

Thalidomide 4'-oxyacetamide-alkylC8-amine

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Functionalized cereblon ligand for PROTAC research and development; incorporates an E3 ligase ligand plus an alkylC8 linker with terminal amine ready for conjugation to a target protein ligand. Part of a range of functionalized tool molecules for PROTAC R&D. This product has been recently renamed. The previous name for this product was Thalidomide amide-alkylC8-amine PROTAC is a registered trademark of Arvinas Operations, Inc., and is used under license.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481