


Stephodeline

Chemical Properties

CAS No. : 56596-12-4
 Formula: C₂₁H₂₇N₅O₅
 Molecular Weight: 373.44
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Stephodeline is a natural product for research related to life sciences. The catalog number is TN6261 and the CAS number is 56596-12-4.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6778 mL	13.389 mL	26.7781 mL
5 mM	0.5356 mL	2.6778 mL	5.3556 mL
10 mM	0.2678 mL	1.3389 mL	2.6778 mL
50 mM	0.0536 mL	0.2678 mL	0.5356 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481