

Compound 0449-0077

Chemical Properties


CAS No. :

Formula: C₂₅H₃₃N₂O₆

Molecular Weight: 443.54

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound 0449-0077 is a useful organic compound for research related to life sciences and the catalog number is T131668.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2546 mL	11.2729 mL	22.5459 mL
5 mM	0.4509 mL	2.2546 mL	4.5092 mL
10 mM	0.2255 mL	1.1273 mL	2.2546 mL
50 mM	0.0451 mL	0.2255 mL	0.4509 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481