

3-[(3-[[3-hydroxy-2-(3-hydroxydecanamido)propanoyl]oxy]decanoyl)oxy]octanoic acid

Chemical Properties

CAS No. :

Formula: C31H57NO9

Molecular Weight: 587.8

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	3-[(3-[[3-hydroxy-2-(3-hydroxydecanamido)propanoyl]oxy]decanoyl)oxy]octanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7013 mL	8.5063 mL	17.0126 mL
5 mM	0.3403 mL	1.7013 mL	3.4025 mL
10 mM	0.1701 mL	0.8506 mL	1.7013 mL
50 mM	0.034 mL	0.1701 mL	0.3403 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481