


## Talastine

## Chemical Properties

CAS No. :	16188-61-7
Formula:	C <sub>19</sub> H <sub>21</sub> N <sub>3</sub> O
Molecular Weight:	307.39
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Talastine is an antihistamine compound that can induce an allergic rash.
Targets(IC50)	Histamine Receptor

## Solubility Information

Solubility	DMSO: 55 mg/mL (178.93 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2532 mL	16.266 mL	32.532 mL
5 mM	0.6506 mL	3.2532 mL	6.5064 mL
10 mM	0.3253 mL	1.6266 mL	3.2532 mL
50 mM	0.0651 mL	0.3253 mL	0.6506 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Richter G, ET AL. Talastine (Ahanon) as a cause of allergic drug exanthema. Dermatol Monatsschr. 1990;176(2-3): 111-113.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481