

## Doronanine

## Chemical Properties


CAS No. :

Formula: C<sub>18</sub>H<sub>25</sub>NO<sub>5</sub>

Molecular Weight: 335.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Doronanine is a useful organic compound for research related to life sciences and the catalog number is T124290.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9815 mL	14.9076 mL	29.8151 mL
5 mM	0.5963 mL	2.9815 mL	5.963 mL
10 mM	0.2982 mL	1.4908 mL	2.9815 mL
50 mM	0.0596 mL	0.2982 mL	0.5963 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481