


## copyranosyl ester

## Chemical Properties

CAS No. : 81263-98-1  
 Formula: C<sub>26</sub>H<sub>38</sub>O<sub>10</sub>  
 Molecular Weight: 510.58  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	copyranosyl ester is a natural product for research related to life sciences. The catalog number is TN3952 and the CAS number is 81263-98-1.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9586 mL	9.7928 mL	19.5856 mL
5 mM	0.3917 mL	1.9586 mL	3.9171 mL
10 mM	0.1959 mL	0.9793 mL	1.9586 mL
50 mM	0.0392 mL	0.1959 mL	0.3917 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481