


KY1220

Chemical Properties

CAS No. : 292168-79-7
Formula: C₁₄H₁₀N₄O₃S
Molecular Weight: 314.32
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	KY1220 destabilizes both β -catenin and Ras by targeting the Wnt/ β -catenin pathway. It has an IC ₅₀ of 2.1 μ M in HEK293 reporter cells.
Targets(IC ₅₀)	Wnt/beta-catenin
In vitro	KY1220 accelerates the degradation rates of both β -catenin and Ras in SW480 cell lines. KY1220 dose-dependently reduces Wnt3a-CM-induced TOPflash reporter activation and mRNA expression of Wnt target genes CCND1 and MYC in HEK293 cells. Ras destabilization by KY1220 consequently inhibits the activities of both ERK and Akt, which are downstream effectors of Ras in SW480 cells harboring a KRAS mutation. In HEK293 cells, both β -catenin and pan-Ras protein levels are similarly reduced in a dose-dependent manner after treatment with KY1220, whereas the mRNA levels of CTNNB1 (which encodes β -catenin), NRAS, KRAS, and HRAS remain unchanged. K-Ras, which has a critical role in the progression of CRCs, is also destabilized by KY1220 via polyubiquitin-dependent proteasomal degradation. The proliferation and transformation of the HCT15, SW480, D-WT, and D-MT CRC cells are efficiently inhibited after treatment with KY1220.

Solubility Information

Solubility	DMSO: 100 mg/mL (318.15 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1815 mL	15.9074 mL	31.8147 mL
5 mM	0.6363 mL	3.1815 mL	6.3629 mL
10 mM	0.3181 mL	1.5907 mL	3.1815 mL
50 mM	0.0636 mL	0.3181 mL	0.6363 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Cha PH, et al. Small-molecule binding of the axin RGS domain promotes β -catenin and Ras degradation. Nat Chem Biol. 2016 Aug;12(8):593-600.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481