


AMPCP

Chemical Properties

CAS No. :	104835-70-3
Formula:	C ₁₁ H ₁₅ N ₅ Na ₂ O ₉ P ₂
Molecular Weight:	469.19
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	AMPCP is an Ecto-5'-nucleotidase (CD73) inhibitor.
Targets(IC50)	Others,CD73

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1313 mL	10.6567 mL	21.3133 mL
5 mM	0.4263 mL	2.1313 mL	4.2627 mL
10 mM	0.2131 mL	1.0657 mL	2.1313 mL
50 mM	0.0426 mL	0.2131 mL	0.4263 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481