

## Compound T130683

## Chemical Properties


CAS No. :

Formula: C45H72O16

Molecular Weight: 869.06

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2S,3R,4R,5R,6S)-2-(((2R,3S,4S,5R,6R)-4-hydroxy-2-(hydroxymethyl)-6-((4'S,7'S,9'S,13'R,16'S)-5,7',9',13'-tetramethyl-5'-oxaspiro[oxane-2,6'-pentacyclo[10.8.0.0 <sup>2,9</sup> .0 <sup>4,8</sup> .0 <sup>13,18</sup> ]icosan]-18'-eneoxy)-5-(((2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)oxan-3-yl]oxy)-6-methyloxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1507 mL	5.7533 mL	11.5067 mL
5 mM	0.2301 mL	1.1507 mL	2.3013 mL
10 mM	0.1151 mL	0.5753 mL	1.1507 mL
50 mM	0.023 mL	0.1151 mL	0.2301 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481