


SQ 31906

Chemical Properties

CAS No. : 81131-74-0
 Formula: C₂₃H₃₆O₇
 Molecular Weight: 424.53
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|---------------------------------------|
| Description | SQ 31906 is a prevastatin metabolite. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.3555 mL | 11.7777 mL | 23.5555 mL |
| 5 mM | 0.4711 mL | 2.3555 mL | 4.7111 mL |
| 10 mM | 0.2356 mL | 1.1778 mL | 2.3555 mL |
| 50 mM | 0.0471 mL | 0.2356 mL | 0.4711 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Mulvana D, Jemal M, Pulver SC. Quantitative determination of pravastatin and its biotransformation products in human serum by turbo ion spray LC/MS/MS. J Pharm Biomed Anal. 2000 Oct;23(5):851-66. PubMed PMID: 11022911.

Jemal M, Xia YQ, Whigan DB. The use of high-flow high performance liquid chromatography coupled with positive and negative ion electrospray tandem mass spectrometry for quantitative bioanalysis via direct injection of the plasma/serum samples. Rapid Commun Mass Spectrom. 1998;12(19):1389-99. Erratum in: Rapid Commun Mass Spectrom 1998;12(23):1978. Yuan-Qing [corrected to Xia YQ]. PubMed PMID: 9773524.

Gehr TW, Sica DA, Slugg PH, Hammett JL, Raymond R, Ford NF. The pharmacokinetics of pravastatin in patients on chronic hemodialysis. Eur J Clin Pharmacol. 1997;53(2):117-21. PubMed PMID: 9403282.

Halstenson CE, Triscari J, DeVault A, Shapiro B, Keane W, Pan H. Single-dose pharmacokinetics of pravastatin and metabolites in patients with renal impairment. J Clin Pharmacol. 1992 Feb;32(2):124-32. PubMed PMID: 1613121.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481