


## 3,5-Di-O-galloylquinic acid

## Chemical Properties

CAS No. : 86588-90-1  
 Formula: C<sub>21</sub>H<sub>20</sub>O<sub>14</sub>  
 Molecular Weight: 496.377  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3,5-Di-O-galloylquinic acid is a useful organic compound for research related to life sciences. The catalog number is T124062 and the CAS number is 86588-90-1.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0146 mL	10.0729 mL	20.1459 mL
5 mM	0.4029 mL	2.0146 mL	4.0292 mL
10 mM	0.2015 mL	1.0073 mL	2.0146 mL
50 mM	0.0403 mL	0.2015 mL	0.4029 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481