


GSK-598809 L-tartrate sesquihydrate

Chemical Properties

CAS No. : 1008521-09-2
 Formula: C₂₆H₂₉F₄N₅O₇S.3/2 H₂O
 Molecular Weight: 658.62
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	GSK-598809 L-tartrate sesquihydrate is an acute D3 antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5183 mL	7.5916 mL	15.1833 mL
5 mM	0.3037 mL	1.5183 mL	3.0367 mL
10 mM	0.1518 mL	0.7592 mL	1.5183 mL
50 mM	0.0304 mL	0.1518 mL	0.3037 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481