


3-(2-Naphthalenyl)-1-(3-pyridinyl)-2-propen-1-one

Chemical Properties

CAS No. :	1031063-36-1
Formula:	C ₁₈ H ₁₃ NO
Molecular Weight:	259.3
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	DMU2105 is a specific CYP1B1 inhibitor(CYP1B1 and CYP1A1 with IC50 of 10 nM and 742 nM, respectively)
Targets(IC50)	Cytochromes P450

Solubility Information

Solubility	DMSO: 25 mg/mL (96.41 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+90% Corn Oil: 2 mg/mL (7.71 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8565 mL	19.2827 mL	38.5654 mL
5 mM	0.7713 mL	3.8565 mL	7.7131 mL
10 mM	0.3857 mL	1.9283 mL	3.8565 mL
50 mM	0.0771 mL	0.3857 mL	0.7713 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Horley N J , Beresford K J M , Chawla T , et al. Discovery and characterization of novel CYP1B1 inhibitors based on heterocyclic chalcones: Overcoming cisplatin resistance in CYP1B1-overexpressing lines[J]. European Journal of Medicinal Chemistry, 2017, 129:159-174.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481