

Fmoc-N-PEG-CH₂COOH (MW 5000)**Chemical Properties**

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fmoc-N-amido-PEG-CH ₂ COOH (MW 5000), a PEG-based PROTAC linker, is used in PROTAC synthesis [1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together; one ligand binds to an E3 ubiquitin ligase, and the other targets a specific protein. By utilizing the intracellular ubiquitin-proteasome system, PROTACs selectively degrade target proteins[1].

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481