

6-S-(2-Cyanoethyl)-2'-deoxy-5'-O-DMTr--6-thioguanosine 3'-CED phosphoramidite

Chemical Properties

CAS No. :


Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	6-S-(2-Cyanoethyl)-2'-deoxy-5'-O-DMTr-6-thioguanosine 3'-CED phosphoramidite [catalog number TNU1665] is an essential organic compound for life sciences research.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481