


## 4-chlorofuro[3,2-c]pyridine

## Chemical Properties

CAS No. :	31270-80-1
Formula:	C7H4ClNO
Molecular Weight:	153.57
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-chlorofuro[3,2-c]pyridine is a molecular block that can be used as a starting material for the synthesis of antibiotics, antitumor compounds and other bioactive substances.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.5117 mL	32.5584 mL	65.1169 mL
5 mM	1.3023 mL	6.5117 mL	13.0234 mL
10 mM	0.6512 mL	3.2558 mL	6.5117 mL
50 mM	0.1302 mL	0.6512 mL	1.3023 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481